


Analysenzertifikat

Certificate of Analysis



virion\serion

SERION ELISA *control* PARVOVIRUS B19 IgM

Best.-Nr. / Order No.	BR122M
Chargen-Nr. / Lot No.	SIZ.BE
Verfallsdatum / Expiry Date	2019-09
OD (405 nm) Zielwert / Target Value	1,160
OD (405 nm) Zielbereich / Target Range	$\geq 0,580$ und / and $\leq 2,320$
Zielwert (Titer) / Target Value (Titre)	56,8 U/ml
Zielbereich (Titer) / Target Range (Titre)	$\geq 28,4$ U/ml und / and $\leq 113,6$ U/ml
Rekonstitution / Reconstitution	0,15 ml Aqua dest.
Serumverdünnung / Serum Dilution	1 + 100 in DILBS3
Bitte RF-Absorbens verwenden! Please use Rf-absorbent!	Bitte Gebrauchsanweisung zu RF- Absorbens Z 4 / Z 200 beachten! Please refer to the instruction manual of Rf-absorbent Z 4 / Z 200!
Versions-Nr. / Version No.	2010 / 03-1

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de