



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION ELISA *control* RESPIRATORY SYNCYTIAL VIRUS (RSV) IgG

Best.-Nr. / Order No.	BR113G
Chargen-Nr. / Lot No.	SHA.CC
Verfallsdatum / Expiry Date	2020-08
OD (405 nm) Zielwert / Target Value	0,813
OD (405 nm) Zielbereich / Target Range	$\geq 0,406$ und / and $\leq 1,626$
Zielwert (Titer) / Target Value (Titre)	24,4 U/ml
Zielbereich (Titer) / Target Range (Titre)	$\geq 12,2$ U/ml und / and $\leq 48,8$ U/ml
Rekonstitution / Reconstitution	0,10 ml Aqua dest.
Serumverdünnung / Serum Dilution	1 + 2000
Versions-Nr. / Version No.	2010 / 09-1

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de