


# Analysenzertifikat

## Certificate of Analysis



virion\sersion

### SERION ELISA *control* VARICELLA-ZOSTER VIRUS IgM

Best.-Nr. / Order No.	BR104M
Chargen-Nr. / Lot No.	SFA.FS
Verfallsdatum / Expiry Date	2020-06
OD (405 nm) Zielwert / Target Value	1,316
OD (405 nm) Zielbereich / Target Range	$\geq 0,658$ und / and $\leq 2,631$
Zielwert (Titer) / Target Value (Titre)	22,9 U/ml
Zielbereich (Titer) / Target Range (Titre)	$\geq 11,5$ U/ml und / and $\leq 45,8$ U/ml
Rekonstitution / Reconstitution	0,1 ml Aqua dest.
Serumverdünnung / Serum Dilution	1 + 100
Bitte RF-Absorbens verwenden! Please use Rf-absorbent!	Bitte Gebrauchsanweisung zu RF- Absorbens T 4 / T 200 beachten! Please refer to the instruction manual of Rf-absorbent T 4 / T 200!
Versions-Nr. / Version No.	2010 / 10-1

Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)