

# Product Data Sheet



## Bulk antigen Echinococcus granulosus

virion\serion

<b>Order no.</b>	BA107VS
<b>Lot no.</b>	Lot specific
<b>Agent</b>	Echinococcus granulosus
<b>Culture system</b>	Material derived from bovine hydatid cyst.
<b>Preparation</b>	Bovine hydatid cyst fluid, sterile filtered.
<b>Inactivation</b>	Whole infective organisms removed by filtration. Inactivity is confirmed by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.
<b>Preservation</b>	None
<b>Protein concentration</b>	Lot specific, specified in mg/mL. Determined by Bradford assay (Biorad) using BSA as standard.
<b>Analysis</b>	By ELISA
<b>Storage</b>	-20°C to -80°C. Avoid repeated freeze-thaw cycles.
<b>Date of manufacture</b>	Lot specific
<b>Expiry date</b>	Lot specific

Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)

FB-IV-062 A