


Product Data Sheet



virion\serion

Bulk antigen Brucella abortus

Order no.	BA116VS
Lot no.	Lot specific
Agent	Brucella abortus (Strain S99 - Weybridge)
Culture system	Incubated hen's eggs
Preparation	Lysate in glycine buffer (pH 9.5) with 2% (w/v) N-octylglucoside
Inactivation	This product has been treated with heat (1 hour, 65°C), phenol (0.5 % (w/w)) and detergent. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.
Preservation	None
Protein concentration	Lot specific, specified in mg/mL. Determined by Bradford assay (Biorad) using BSA as standard.
Analysis	By ELISA
Storage	-20°C to -80°C. Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific
Expiry date	10 years from manufacture

Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-062 A