

Product Data Sheet



Bulk antigen Helicobacter pylori

virion\sersion

Order no.	BA118VS
Lot no.	Lot specific
Agent	Helicobacter pylori
Culture system	Solid medium
Preparation	Washed harvested organisms are solubilised and partially purified. Antigen is presented in a 1% non-ionic detergent solution.
Inactivation	Detergent treatment. Inactivity is confirmed by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.
Preservation	None
Protein concentration	Lot specific, specified in mg/mL. Determined by Bradford assay (Biorad) using BSA as standard.
Analysis	By ELISA
Storage	-20°C to -80°C. Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific
Expiry date	10 years from manufacture

Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-062 A