

Product Data Sheet



Bulk antigen Influenza B Virus

virion\sersion

Order no.	BA1232VS
Lot no.	Lot specific
Agent	Influenza B virus - Strain Harbin
Culture system	Incubated hen's eggs
Preparation	Surface antigens of the viruses are removed by deoxycholate treatment and the remaining core particles are concentrated by centrifugation. Some surface proteins may be present.
Inactivation	Deoxycholate treatment. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.
Preservation	None
Protein concentration	Lot specific, specified in mg/mL. Determined by Bradford assay (Biorad) using BSA as standard.
Analysis	By ELISA
Storage	-20°C to -80°C. Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific
Expiry date	10 years from manufacture

Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-062 A