

# Product Data Sheet



## Bulk antigen Aspergillus fumigatus

virion\sersion

<b>Order no.</b>	BA132VS
<b>Lot no.</b>	Lot specific
<b>Agent</b>	Aspergillus fumigatus
<b>Culture system</b>	Liquid culture
<b>Preparation</b>	Disrupted somatic material clarified by centrifugation in phosphate buffered saline.
<b>Inactivation</b>	Viable organisms are removed by centrifugation. Inactivity is confirmed by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.
<b>Preservation</b>	None
<b>Protein concentration</b>	Lot specific, specified in mg/mL. Determined by Bradford assay (Biorad) using BSA as standard.
<b>Analysis</b>	By ELISA
<b>Storage</b>	-20°C to -80°C. Avoid repeated freeze-thaw cycles.
<b>Date of manufacture</b>	Lot specific
<b>Expiry date</b>	10 years from manufacture

Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)

FB-IV-062 A