

Product Data Sheet



virion\sersion

Bulk antigen Yersinia enterocolitica O:9 YOP

Order no.	BA13809VS
Lot no.	Lot specific
Agent	Yersinia enterocolitica serovar O:9
Culture system	Liquid culture
Preparation	Yersinia is grown in liquid culture and stimulated to secrete YOPs (Yersinia outer proteins). YOP proteins concentrated from the culture medium are presented in phosphate buffered saline.
Inactivation	Sterile filtration. Inactivity is confirmed by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.
Preservation	None
Protein concentration	Lot specific, specified in mg/mL. Determined by Bradford assay (Biorad) using BSA as standard.
Analysis	By SDS-PAGE and ELISA
Storage	-20°C to -80°C. Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific
Expiry date	10 years from manufacture

Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-062 A