

Product Data Sheet



virion\sersion

Bulk antigen Campylobacter jejuni

Order no.	BA139VS
Lot no.	Lot specific
Agent	Campylobacter jejuni
Culture system	Blood agar
Preparation	Partially purified detergent extract of the membrane fraction. The antigen is presented in water. The main component is an outer membrane protein with approximately 45 kDa in SDS-PAGE.
Inactivation	Disruption and detergent treatment. Inactivity is confirmed by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.
Preservation	None
Protein concentration	Lot specific, specified in mg/mL. Determined by Bradford assay (Biorad) using BSA as standard.
Analysis	By SDS-PAGE and ELISA
Storage	-20°C to -80°C. Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific
Expiry date	10 years from manufacture

Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-062 A