Product Data Sheet – Recombinant Antigen

Name of product	PR3
Order No.	BA502R01
Lot No.	Lot specific
Source	Recombinant protein, expressed in insect cells
Product description	Recombinant human PR3 antigen (Leukocyte proteinase 3, c-ANCA antigen, Wegener autoantigen); Molecular mass: 26 kDa
Affinity tag	Strep-tag
Purification	By chromatography
Presentation	The protein is presented in 50 mM HEPES/NaOH pH 8.0, 150 mM NaCl, 10 % glycerol and 0.0098% ProClin950
Storage	-65 °C or below. Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific
Expiry	10 years from date of manufacture (when stored under recommended conditions)
Manufacturer	Institut Virion\Serion GmbH, Germany

For research or manufacturing use only. Not for use in diagnostic procedures. 2025-12-15

Analysis:

Test	Method	Specification
Protein concentration	Bradford assay (Bio-Rad) with BSA as reference standard	0.5 – 2 mg/mL
Identity / Purity	SDS-PAGE (reducing)	Major bands comparable to Reference Standard
Immunoreactivity	Bead-based Immunoassay IgG	QC Comparison to master lot within specified range

virion\serion

Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19, 97076 Würzburg, Germany
Phone +49 931 3045 0 **Fax** +49 931 3045 100 **Mail** immunologics@virion-serion.de **Web** www.virion-serion.de